


CHEMICAL WASTE COMPATIBILITY TEST

Chemical waste reactions can be complex and sometimes delayed. This test provides a quick check on compatibility, but further judgment is required. Always consult your supervisor or HSEO if unsure.





Safety First

- Wear safety goggles, gloves, and a lab coat
- Conduct the test inside a fume cupboard
- Do NOT conduct the test alone




Prepare the Sample

- Homogenize the chemical waste inside the target container
- Withdraw 50 mL into a beaker with a minimum capacity of 250 mL
- Measure & record the temperature



Mix with New Waste

- Slowly add a portion of the new waste in proportion to actual disposal volumes (e.g., 10 mL for a 1:5 ratio)
- STOP IMMEDIATELY if a rigorous reaction occurs!



Observe & Measure

- After mixing, re-measure the temperature
- If **bubbling**, **gassing**, or a temperature rise of ≥10°C occurs → INCOMPATIBLE!
- Store separately & label with chemical composition, date, and producer contact information

Final Assessment

- Check for colour change or precipitation alone. Compatible if NO other adverse reactions occur within 5 minutes
- If incompatible, pre-treat waste before disposal or use the **Labpack Waste Scheme**