# Safe Operating Procedure of using Carbon Tetrachloride for Animal Research

NOTE: You must read this entire document and both you and the Principal Investigator must sign it before commencing any work.

Project Number (e.g. Tick@Lab):{ University SOP}	
Principal Investigator/Supervisor:{ University SOP}	
Room and Building where SOP is used: Wet labs, LAF labs	

# Summary of How the Material/Equipment/Process will be Used

Carbon tetrachloride is commonly used in experimental animals to induce chronic liver injury.

The processes associated with the use of carbon tetrachloride include:

- 1) Preparation of Carbon tetrachloride solution
- 2) Transport of Carbon tetrachloride from research lab to LAF for administration
- 3) Administration of Carbon tetrachloride into animals
- 4) Housing and handling of Carbon tetrachloride administered animals
- 5) Cage changing and washing of Carbon tetrachloride contaminated cages
- 6) Disposal of Carbon tetrachloride contaminated waste

### **Potential Hazards**

Human Toxicity or Carcinogenic Exposure (Process 1-6)

- Toxic if swallowed, inhaled or in contact with skin (hazard statement H301, H331 and H311)
- May cause cancer (hazard statement H351)
- Damages the specific organ through prolonged or repeated exposure (hazard statement H372)

Needle Stick Injury due to Injection (Process 3)

Animal Bite (Process 3-5)

Carbon tetrachloride SDS weblink:

https://www.sigmaaldrich.com/HK/en/sds/sial/289116?userType=anonymous

#### **Safety Installations**

# **Engineering Control Measures**

- Certified fume hood (FH) for Carbon tetrachloride preparation (e.g. open, weight, reconstitute and dilute)
- Certified Class II Biosafety cabinet (BSC) for injection, animal handling and cage changing
- Secondary container for transport
- Luer lock syringes for administration
- Animal restraint devices for administration (if applicable)
- Forceps for animal handling
- Animal cage changing station for cage changing
- Individually ventilated cage (IVC) and air handling unit for housing
- Cage liners for housing
- Cage washing machines for washing
- Carbon tetrachloride waste bin for disposal of contaminated waste
- Carbon tetrachloride sharp box for disposal of contaminated syringes
- Eye wash station for emergency use
- Availability of Carbon tetrachloride Safety Data Sheet (SDS) in handling areas

### **Work Practices**

### **LAF Notification**

• Inform LAF prior to using Carbon tetrachloride in LAF for better arrangement

## Medical Surveillance Program

• Enroll in medical surveillance program and indicate "toxic chemicals injected"

# **Training**

- Complete Carbon tetrachloride-specific safety training provided by the PI
- Complete HSEO MC03 Chemical I, MC07 Chemical II, MC06 Biological Safety Training
- Complete general animal user training and animal handling training provided by LAF

### SOP

• Establish and follow in-house SOPs for specific process (e.g. Carbon tetrachloride administration, waste collection, cage washing, etc.)

### **Labelling and Warning Sign**

- Complete LAF blue cage card "Health Hazard Card" with administration information and place in each Carbon tetrachloride administered animal cage until 72 hours post-administration.
- Affix a Carbon tetrachloride warning sign with hazard information on each of the Carbon tetrachloride administered animal holding room
- Label all Carbon tetrachloride tubes, containers, waste bin, waste bag with appropriate hazard warnings and chemical identification.
- Post a "Carbon Tetrachloride Hazard" sign on the FH and BSC which is used for Carbon tetrachloride preparation and administration

#### **Proper Storage**

- Store Carbon tetrachloride in a lockable secondary container, and in a well-ventilated place
- Use zipper bags and secondary containers during transport to avoid spillage

# Good Hygiene Practice

- Minimize exposure by implementing control measures and wear adequate PPE
- Clean work surfaces with detergent and water, followed by 10% sodium hypochlorite and rinse thoroughly
- Wash hands immediately after handling Carbon tetrachloride and Carbon tetrachloride-contaminated waste

# **Specific Experimental Procedures**

Preparation, Transportation and Administration of Carbon Tetrachloride (process 1-3)

- Purchase only the necessary quantity of Carbon tetrachloride
- Wear adequate PPE when handling Carbon tetrachloride
- Inspect the outer package of Carbon tetrachloride to make sure it is intact
- Open, weigh and prepare Carbon tetrachloride solution inside a certified FH
- A plastic-backed absorbent pad should be placed under the work surface during preparation and administration process to avoid contamination
- Place a waste bag on the work surface for the collection of contaminated wastes
- Aliquot the Carbon tetrachloride solution into leak-proof, screw cap tubes and place inside a secondary container labelled with chemical name and hazard warning labels. Store the container in a fridge.
- During transportation, Carbon tetrachloride solution tubes must be stored in a zipper bag inside a secondary container (only bring the required amount to LAF)
- Conduct the administration inside a certified BSC
- Do not recap the syringes. Dispose of syringes used into Carbon tetrachloride sharp box inside BSC
- Dispose of all contaminated wastes (e.g. tubes, tips, absorbent pad, syringes, etc.) into a waste bag inside FH/BSC and then discard into designated Carbon tetrachloride chemical waste bin
- Clean work surfaces with detergent and water, followed by 10% sodium hypochlorite and rinse thoroughly.

Housing and Handling of Carbon Tetrachloride-Administered Animals (process 4)

- Complete LAF blue cage card "Health Hazard Card" with administration information and place in each Carbon tetrachloride administered animal cage
- Conduct cage changes within a BSC or Cage Changing Station (CCS)
- Each Carbon tetrachloride-administered animal should return to a cage with the cage liner and "Health Hazard Card", then housed in a designated animal holding room

Cage Changing and Washing of Carbon Tetrachloride-Contaminated Cages (process 5)

- After 72 hours of last administration, the contaminated cage/water/feed can be changed inside a BSC or CCS
- Dispose all contaminated disposable waste into the cage liner inside BSC/CCS, and bag the cage liners into a chemical waste bag (refer to waste disposal for more details)
- After disposing wastes, the "Health Hazard Card" can be removed and considered free of Carbon tetrachloride
- Bag the non-disposable cage accessories' (e.g. lid, cage base, wire bar, etc.) in a bag for LAF collection
- LAF staff collect the bagged cage accessories to cage washing area for normal washing process

## Personal Protective Equipment

Skin Protection, Eye Protection, Face Protection and Respiratory Protection

- Wear double-layer Neoprene gloves. One glove should be placed under the coat / coveralls and one
  over.
- Lab coat or coveralls
- Sleeve covers if wrist is exposed
- Safety glasses or chemical goggles or face shield
- N95 mask or PAPR if work is not conducted in BSC or fume hood

# **Waste Disposal**

## Carbon Tetrachloride-Contaminated Waste Disposal (process 6)

Discard all contaminated wastes during preparation and administration into waste bags.

### Double bags all wastes.

- Tie up the waste bags and place in designated Carbon tetrachloride chemical waste bin for HSEO collection
- The waste bin should be kept closed at all times
- For contaminated sharps waste box, place the sharp box in designated Carbon tetrachloride chemical waste bin for HSEO collection
- For Carbon tetrachloride-administered animal carcasses, **double bags** before placing in the clinical waste fridge for contractor collection

# Bedding and Cage Liner Disposal (process 6)

- · Place the bedding and cage liners in a chemical waste bag, and tie up the bag
- · Place the waste bag in designated Carbon tetrachloride chemical waste bin for HSEO collection

# **Spills and Incidents**

# Minor Spill of Carbon Tetrachloride

- · Absorb the spill with absorbent pad or paper towel with adequate PPE
- Cover, Clean work surfaces with detergent and water, followed by 10% sodium hypochlorite and rinse thoroughly
- Collect the wastes in a plastic bag and dispose in a Carbon tetrachloride-contaminated waste bin

# Major Spill of Carbon Tetrachloride

- · Avoid inhalation and generating dust
- · Remove contaminated clothing with gloved hands, remove gloves and place near the spill
- Remove all ignitable sources. If a fire is visible, DO NOT turn on the emergency ventilation.
- · Evacuate other workers within the laboratory.
- Call Security Control Center Ext 8999 for assistance and do not allow anyone to enter the affected area
- · HSEO will be notified by Security Control Center to clean up the spill

### Bites by Animals Treated with Carbon Tetrachloride

- Wash bite area with running water for 15 minutes
- Seek medical attention. Bring SDS to the clinic / doctor.

# **Emergency Procedures**

# **Exposure**

- If Carbon tetrachloride comes into contact with the eye, face and skin, wash the area with running water for 15 mins using the nearest emergency facilities and seek medical advice immediately.
- · Notify HSEO of any exposures.

# **Emergency Response**

• Notify the Security Control Center by dialing ext. 8999 and provide information on the incident, including the chemical involved, the location, and any injuries.

#### References

- Benson, J. M., & Springer, D. L. (1999). Improved risk estimates for carbon tetrachloride (No. EMSP-54940-). Lovelace Biomedical and Environmental Research Institute, Albuquerque, New Mexico; Pacific Northwest National Lab.(PNNL), Richland, WA (United States).
- Cho KJ, Jang JJ. Effects of carbon tetrachloride, ethanol and acetaldehyde on diethylnitrosamine-induced hepatocarcinogenesis in rats. Cancer Lett. 1993 Jun 15;70(1-2):33-9. doi: 10.1016/0304-3835(93)90071-g. PMID: 8330298.
- National Center for Biotechnology Information (2025). PubChem Compound Summary for CID 5943, Carbon Tetrachloride. Retrieved November 5, 2025 from <a href="https://pubchem.ncbi.nlm.nih.gov/compound/Carbon-Tetrachloride">https://pubchem.ncbi.nlm.nih.gov/compound/Carbon-Tetrachloride</a>.

Prepared by:	Date:November 12 <sup>th</sup> , 2025	_November 12 <sup>th</sup> , 2025
Reviewed by:	Date:	Date: